



### 1080P Wi-Fi Digital Microscope for iOS/Android/PC, 10x-230x Magnification

Model No: UM018G

- \* 1080P Full HD Wi-Fi to iOS, Android device
- \* metal stand with sample clips and object focusing
- \* Video & Photo capture
- \* Rechargeable Li-ion Battery, SD card slot and 16GB SD card
- \* Also can be used as wired by USB to computer
- \* Measurement available when as USB microscope

\* **Magnifications:** 10x-230x

\* **Focal Distance:** 0mm to 50mm

\* **Camera Type:** high quality CMOS

\* **Video Resolution:** 1080P (1920x1080), 720P (1280x720), VGA (640x480)

\* **Video resolution:** 30FPS

\* **Picture Resolution:** 3M, 5M, 8M, 12M \***Resolution:** 3M

\* **Picture Format:** JPG

\* **Wireless:** Wi-Fi

\* **Transmission Range:** 20M (open sight)

\* **Transmission Rate:** 150Mbps,

\* **Power supply:** 1) Built-in rechargeable lithium battery; 2) DC 5V/1A (mini USB)

\* **Battery life:** around 100minutes

\* **Battery charging time:** around 90 minutes

\* **Measurement:** available on PC (by USB connection)

\* **Operation Temperature:** 0-55 degree

\* **Light source:** 6 LEDs (brightness adjustable)

\* **System supported:** iOS (5.1 or later), Android, Windows-XP (SP2), W7/8/10 compatible software